


# DOTZEN GANGLIO PROFILE Ab




ZenTech

Immunodot assay for the determination of Ig(G+M) and/or IgG & IgM autoantibodies specific to **gangliosides** in human serum

Dot Blot immunoassay (24 determinations)

DOTZEN Ganglio Profile Ab Code: L-ZA-24



Getting better  
**NEUROPATHIES**  
diagnosis using the  
10 most important  
epitopes

DOTZEN Ganglio Profile Ab is a useful help to :

- **Confirm** & establish a diagnosis of neuropathy
- **Discriminate** between neuropathies and other etiologies
- **Identify** and classify the anti-gangliosides antibodies

GANGLIO PROFILE	GANGLIO PROFILE
Sulfatides	●
GQ1b	●
GT1b	●
GT1a	●
GD3	●
GD1b	●
GD1a	●
GM3	●
GM2	●
GM1	●
Funct Ctl	●

- Easy and efficient DOT technique
- Allows testing of the most important antigens simultaneously
  - 10 antigens spotted on PVDF membrane
  - Large strip allowing easy reading
  - Stability of coloration
- First step : screening combining the detection of IgG & IgM isotypes
- Second step : IgG and IgM identification (enzyme conjugates in the same kit)
- 1/505 final sera dilution

Liège Science Park  
Avenue du Pré Aily ,10  
4031 Angleur – BELGIUM

Tel: +32 (0)4 3614232  
Fax: +32 (0)4 3670063

info@zentech.be  
[www.zentech.be](http://www.zentech.be)

## What are gangliosides ?

Gangliosides are glycosphingolipids with sialic acid. Gangliosides differ from one other according to the number and position of the hexoses and the residues of sialic acid . Molecular targets are thus numerous and varied, giving rise to the existence of numerous antibodies of different specificities. The nature, quantity and distribution of the diverse gangliosides in the nervous system vary. The motor nerves, for instance, are rich in GM1 and GD1b gangliosides while the sensory nerves mainly contain polysialogangliosides. **There is a good correlation between the type of gangliosides recognized by the antibodies and the nature of the neurological impairment.**


## Search for anti-ganglioside antibodies

Due to the diversity of anti-ganglioside antibodies, it is important to test serums on a large range of gangliosides and to determine which IgG and IgM isotypes they belong to. Immobilization of gangliosides on the membranes of polyvinylidene difluoride (PVDF) used in the DOT technique is easy and efficient. **Immunodot is the preferred technique for the search for anti-gangliosides because it allows testing of serums with different gangliosides simultaneously and establishing a complete profile of antibodies, all in the same step.** Prolonged incubations with gangliosides are necessary given the low affinity of anti-ganglioside antibodies. Anti-gangliosides, on the other hand, are antibodies which react best to cold, which implies that reactions must be carried out at the temperature of 4°C.

## Principle of the DOTZEN test

The **DOTZEN Ganglioside Profile Ab** is an in vitro qualitative determination kit for IgG or IgM anti-ganglioside antibodies in human serum for 10 types of gangliosides (see table) .

**In this way this immuno-dot has proved to be the most sensitive and specific method.**



## Interpretation of positive results\*

ANTIBODIES to	ISOTYPES	ASSOCIATED DISEASES
<b>GM1</b> (monospecific or + GD1b)	IgM	Multifocal Motor Neuropathy (MMN) - Pure Motor, Assymetric, Distal, Upper Limbs, Areflexia, Fasciculations - Motor Conduction Blocks frequent
<b>GM1</b>	IgM	Partial Epilepsia
<b>GM1</b>	IgG	Acute Motor Axonal Neuropathy (AMAN) Acute paralysis, severe, inexcitability of motor nerves
<b>GM2</b> (monospecific)	IgM	Chronic Inflammatory Demyelinating Polyneuropathy
<b>GM2</b>	IgG	Acute Inflammatory Demyelinating Polyneuropathy - Sensori-Motor, Flacid Paralysis, Areflexia - frequent CMV Infections
<b>GM3</b>	IgM	Chronic Inflammatory, Sensori-motor Polyneuropathy - Proximal > Distal, Symmetric, Areflexia
<b>GD1a</b> (monospecific)	IgM	Chronic Motor Neuropathy - Distal, Asymmetric, Areflexia
<b>GD1a</b>	IgG	Subacute Severe Axonal Motor Neuropathy
<b>GD1b</b> (monospecific)	IgM	Relapsing Sensory Ataxic Neuropathy - Chronic, Distal, Areflexia, Ataxia
<b>GD1b</b>	IgG	Acute Sensory Ataxic Neuropathy - Distal Paresthesia, Areflexia, Ataxia
<b>GD3</b>	IgM	Chronic Ataxic Neuropathy
<b>GT1a</b> (monospecific)	IgG	Pharyngo-Cervical-Brachial Palsy - Facial Palsy, Dysphagia, Respiratory failure
<b>GQ1b</b> (cross-reactivity with GT1a)	IgG	Miller-Fisher Syndrome (MFS) - Extraocular muscle Paralysis, Ataxia, Areflexia Bickerstaff's Brain Stem Encephalitis - Coma, Ophtalmoparesis, Ataxia
<b>GQ1b</b> (+ GD1b, GD3, GT1b)	IgM	CANOMAD* - Chronic Ataxic Neuropathy - Monoclonal IgM Gammopathy - Cold Agglutinins
<b>Sulfatide</b>	IgM	Sensory or Sensori-motor Polyneuropathy - Distal, symmetric, Ataxia , Areflexia

\* R-L. Humbel, Laboratoire Luxembourgeois d'Immunopathologie , Luxembourg.

\*\*CANOMAD =Chronic Ataxic Neuropathy Ophtalmo- plegia M-protein Agglutinin Demyelination

